otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

Batch no. 607-03_25315.docx

CE2777

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GB • User instructions**Brand name / type no.****Nilex Ultra PF**

607 070	Size 7
607 080	Size 8
607 090	Size 9
607 100	Size 10
607 110	Size 11

Description

Thin nitrile disposable glove with beaded cuff. The glove is powder free and ambidextrous.

Generally

Before use, it should be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work.

The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

Category

The glove is certified in category III in compliance with Medical Device Regulation 2017/745 and PPE Regulation (EU) 2016/425 on safety requirements for personal protective equipment.

The glove has been tested in accordance with the standards EN 420:2003+A1:2009 (general requirements), EN 455 (medical use) and EN ISO 374-1:2016, EN 374-4:2013, EN ISO 374-5:2016 (chemicals and micro organisms).

Notified body responsible for EU Type Examination and Module C2 ongoing conformity assessment: af SATRA Technology Europe Ltd, Bracetown Business Park, Clonree, Dublin, D15 YN2P, Ireland (Notified Body No. 2777).

EN ISO 374-1:2016/Type B

KPT

EN ISO 374-5:2016

VIRUS

EN 374 Test Data / Type B

EN ISO 374-1:2016

	Chemical	Time*	Level	Degradation %
A	Methanol	N/A	N/A	N/A
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A
I	Ethyl acetate	N/A	N/A	N/A

EN 374-4:2013

J	N-heptane	N/A	N/A	N/A
K	Sodium hydroxide 40%	>480	6	-42.9
L	Sulphuric acid 96%	N/A	N/A	N/A
M	Nitric acid 65 %	<10	0	97.6
N	Acetic acid 99 %	<10	0	93.9
O	Ammonium hydroxide 25%	>10	1	-52.0
P	Hydrogen peroxide 30%	>30	2	22.8
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyde 37 %	>60	3	5.0
	**4% Chlorhexidine Digluconate	>480	6	19.0
	10-13% Sodium Hypochlorite	>480	6	14.7
	50% Sulphuric acid	>480	6	-20.5
	10% Acetic acid	>120	4	66.7
	5% Ethidium Bromide	>480	6	3.4
	50% Glutaraldehyde	>480	6	27.4
	0.1% Phenol	>480	6	33.8
	1.5% Methanol in water	>480	6	21.9
	70% Isopropanol	<10	0	62.2
	35% Ethanol	<10	0	38.8
	3% Povidone-iodine	>480	6	33.7
	10% Sodium Percarbonate	>480	6	15.4

*)Breakthrough time (minutes)

**)Permeation rate: 7 µg/cm²/min

Performance Level	1	2	3	4	5	6
Minimum break-through time (mins)*	10	30	60	120	240	480

*) Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to other factors influencing the performance such as temperature, abrasion, degradation etc.

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EN ISO 374-5:2016

Resistance to Bacteria and Fungi = Pass
Resistance to Virus = Pass

Migration test
(EN 1186-1:2002)**Application / qualities**

Thin and flexible seamless glove with great dexterity. The glove is waterproof and suitable for handling of many different chemicals with short contact time. If the glove has been in contact with substances it must be discarded, when the breakthrough time is reached.

otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

Batch no. 607-03_25315.docx

Primarily used in food production, mounting, product control, laboratories, chemical industry and for general cleaning.

Treatment / storage

The glove is for single use and should not be cleaned or re-used. Best stored in dark and cool surroundings in the original packaging.

Warning

- This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective glove may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- Before usage, inspect the gloves for any defect or imperfections.
- The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimens.

Packaging100 gloves per dispenser.
10 dispensers in a carton made of recyclable cardboard.**SE • Bruksanvisning****Varumärke / art. nr.****Nilex Ultra PF**

607 070	Storlek 7
607 080	Storlek 8
607 090	Storlek 9
607 100	Storlek 10
607 110	Storlek 11

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EN 374-

4:2013

	Chemical	Time*	Level	Degradation %
A	Methanol	N/A	N/A	N/A
B	Acetone	N/A	N/A	N/A
C	Acetonitrile	N/A	N/A	N/A
D	Dichloromethane	N/A	N/A	N/A
E	Carbon disulphide	N/A	N/A	N/A
F	Toluene	N/A	N/A	N/A
G	Diethylamine	N/A	N/A	N/A
H	Tetrahydrofurane	N/A	N/A	N/A

otto schachnerPrinsessens Kvarter 2
7000 Fredericia, DenmarkEU Declaration of conformity available for download at
<http://www.os-safetycenter.com>

Batch no. 607-03_25315.docx

CE2777

I	Ethyl acetate	N/A	N/A	N/A
J	N-heptane	N/A	N/A	N/A
K	Sodium hydroxide 40%	>480	6	-42.9
L	Sulphuric acid 96%	N/A	N/A	N/A
M	Nitric acid 65 %	<10	0	97.6
N	Acetic acid 99 %	<10	0	93.9
O	Ammonium hydroxide 25%	>10	1	-52.0
P	Hydrogen peroxide 30%	>30	2	22.8
S	Hydrofluoric acid 40 %	N/A	N/A	N/A
T	Formaldehyde 37 %	>60	3	5.0
	**4% Chlorhexidine Digluconate	>480	6	19.0
	10-13% Sodium Hypochlorite	>480	6	14.7
	50% Sulphuric acid	>480	6	-20.5
	10% Acetic acid	>120	4	66.7
	5% Ethidium Bromide	>480	6	3.4
	50% Glutaraldehyde	>480	6	27.4
	0.1% Phenol	>480	6	33.8
	1.5% Methanol in water	>480	6	21.9
	70% Isopropanol	<10	0	62.2
	35% Ethanol	<10	0	38.8
	3% Povidone-iodine	>480	6	33.7
	10% Sodium Percarbonate	>480	6	15.4
T	Formaldehyde 37 %	>44	2	33.7

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Description

Green nitrile glove with grip pattern in the palm and on the fingers. Unsupported glove with cotton flock lining inside.

Generally

New and used gloves should be thoroughly inspected before being worn to ensure no damage is present. Before use, it should also be tested/ensured that the glove has the appropriate size to achieve the best possible comfort and safety at work. The service life cannot be determined and depends on the scope of application and the extent to which the user makes sure that the glove is suitable for the intended use.

Category

The glove is certified in category III in accordance with the European PPE Regulation EU 2016/425 on safety requirements for personal protective equipment. The glove has been tested in accordance with the standards EN 420:2003+A1:2009 (general requirements), EN 388:2016+A1:2018 (mechanical) and EN ISO 374-1:2016+A1:2018, EN 374-4:2013, EN ISO 374-5:2016 (chemicals and micro organisms).

CE certified by SATRA Technology Europe Limited, Bracetown Business Park, Clonee, D15 YN2P, Ireland. (notified body number 2777)

Module D assessment by: SGS Fimko Oy, P.O. Box 30 (Särkiniementie 3), 00211 Helsinki, Finland (notified body number 0598)

EN388:2016 EN ISO 374-1:2016/ Type A EN ISO 374-5:2016

4101X AJKLMNOPT

EN 388:2016 Performance Data

Abrasion resistance Level 4 (Max 4)
Blade cut resistance Level 1 (Max 5)
Tear resistance Level 0 (Max 4)
Puncture resistance Level 1 (Max 4)
TDM cut resistance X (Not tested)
Dexterity Level 5

The glove is approved for food handling.

EN ISO 374-1:2016 performance data / Type A :

EN ISO 374-1:2016 Type A EN374-4:2013

	Chemical	Time*	Level	Degradation %
A	Methanol 67-56-1	>30	2	70.1
B	Acetone 67-64-1	N/A	N/A	N/A
C	Acetonitril 75-05-8	N/A	N/A	N/A
D	Dichlormethane 75-09-2	N/A	N/A	N/A
E	Carbon disulphide 75-15-0	N/A	N/A	N/A
F	Toluene 108-88-3	>10	1	81.8
G	Diethylamine 109-89-9	N/A	N/A	N/A
H	Tetrahydrofurane 109-99-9	N/A	N/A	N/A
I	Ethyl acetate 141-78-6	N/A	N/A	N/A
J	N-heptane 142-82-5	>480	6	0.0
K	Sodium hydroxide 40% 1310-73-2	>480	6	4.5
L	Sulphuric acid 96% 7664-93-9	>60	3	61.9
M	Nitric acid 65 % 7697-37-2	>30	2	98.7
N	Acetic acid 99 % 64-19-7	>60	3	91.9
O	Ammonium hydroxide 25% 1336-21-6	>480	6	-5.8
P	Hydrogen peroxide 30 % 7722-84-1	>480	6	-11.7
S	Hydrofluoric acid 40 % 7664-39-3	N/A	N/A	N/A
T	Formaldehyde 37 % 50-00-0	>480	6	-15.6

*Breakthrough time in minutes

otto schachner
Prinsessens Kvarter 2
7000 Fredericia, Denmark
www.os-safetycenter.com (Declaration of Conformity)

Batch no. 900-05_26006

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.

The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm – where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture.

It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

EN 374-4:2013

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

EN ISO 374-5:2016

The penetration resistance has been assessed under laboratory condition and relates only to the tested specimen. This product is not tested for viruses.

Application / qualities

Industry, fishing, agriculture, crafts, food industry and heavy duty cleaning.
Nitrile glove, which provides good protection when handling liquids.

Treatment / storage

Dirty gloves can be cleaned in lukewarm soap water. If the glove has been in contact with chemicals, it is to be discarded once the penetration time limit has been reached.
Chemical cleaning is not recommended.
Must be stored in a dark, dry and cool room in the original packaging. Avoid storage in direct sunlight.

Warning

Gloves are not to be worn when there is a risk of entanglement by moving parts of machines.
This product contains Nitrile Latex and compounding chemicals which may cause allergic reaction in some individuals.

Declaration of conformity is available at
www.os-safetycenter.com

Packaging

1 pair in a PE bag.
12 pairs in a PE bag.
144 pairs in a carton made of recyclable cardboard.

CE 0598

otto schachner
member of ENVA Group

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GB - INSTRUCTIONS FOR USE

Single-use filtering face-piece respirators

IMPORTANT:

Single-use respirators are classified according to one of three categories-FFP1 FFP2, FFP3. Before reading the following information, check the category to which the respirator belongs-this is indicated on the packaging and on the respirator.

If the mask is marked D, it has passed the dolomite clogging test and suitable for use in very dusty environment.

It is the user's responsibility to ensure that the respirator provides the necessary level of protection for the type and concentration of the contaminant(s) in the area where the respirator is intended for use.

DIRECTIONS FOR USE:

- Fit the respirator and check leak-tightness BEFORE entering the contaminated area.
- Wear the respirator for the whole duration of exposure to contaminants.
- Use the respirator in accordance with applicable health and safety regulations.

Discard the respirator and replace with a new one if:

- the respirator is removed while in a contaminated area
- excessive clogging of the respirator causes breathing difficulty or discomfort
- the respirator becomes damaged
- (for respirators protecting against vapours)-the smell of vapours present becomes detectable.

- Leave the contaminated area if dizziness, irritation or other distress occurs.
- Only for single use, No maintenance necessary. Do not store and/or re-use after single use.
- Discard the respirator after single use (one shift, maximum).
- Keep un-used respirators in their closed box and store in a dry non-contaminated area between +5°C and +55°C.

NOTIFIED BODY for CE type examination and for controle of PPE manufactured: Instytut Technologii Tekstylnych CERTEX Sp. z o.o, ul. Gornicza 30/36, 91-765 Lodz, Poland

FITTING INSTRUCTIONS:

1. Cup the respirator in the hand, allowing the head-bands to hang freely.
2. Hold the respirator under chin with the nose-piece facing outwards.
3. Place the lower head-band around the neck below the ears. Holding the respirator against the face with one hand, place the top head-band above the ears, around the crown of the head.
4. Respirators with adjustable head-bands only: while maintaining the position of the respirator, adjust tension by pulling tab of each head-band, (Tension may be decreased by pushing out on the back of the buckle).
5. Mould the nose area to the shape of the face, running the fingertips of both hands from the top of the nose-piece down both sides while pressing inward.

CHECK FOR LEAK-TIGHTNESS (FACIAL FIT) AS FOLLOWS:

- place both hands over the respirator and exhale sharply
- positive pressure should be inside the respirator

- if you feel air escaping around the edges, re-adjust the respirator by tightening the nose-piece and/or increasing tension of the head-bands
- repeat the procedure until satisfactory face-fit has been achieved.

IMPORTANT: In the case of respirators with non-adjustable head-bands-if you feel that the head-bands are too loose, you can increase tension by tying a small knot in the head-band.

LIMITATIONS TO USE:

Do not use these respirators or enter or stay in an area where:

- Oxygen concentration is less than 17% (UK-19%)
- Contaminants or their concentrations are unknown or immediately dangerous to life or health
- Particulate concentrations exceed levels fixed by applicable health and safety regulations or MAC x NPF-whichever is lower, (Nominal protection factors: 4.5 for FFP1, 12.5 for FFP2 and, 50 for FFP3)
NPF: nominal protection factor
MAC: maximum admissible concentration
- Gases and/or vapours are present-except for respirators designated as protecting against vapours, in which case vapour concentrations should not exceed NPF. Respirators protecting against organic vapours should not be used as protection against acid vapours and vice-versa.

WARNING:

Do not use for fire-fighting. These respirators do not supply oxygen. Do not use in oxygen-deficient atmospheres-e.g. tanks or other poorly-ventilated areas (see 'Limitations to use'). Do not use in explosive atmospheres. For use only by trained and qualified personnel. Requirements for leak-tightness are unlikely to be achieved if the respirator is worn against a beard or facial stubble.

When transporting this product, retain in the original packaging and keep away from mechanical and chemical hazards.

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